

兔抗 NGF 多克隆抗体

中文名称：兔抗 NGF 多克隆抗体

英文名称：Anti-NGF rabbit polyclonal antibody

别 名：NGFB; HSAN5; Beta-NGF

抗 原：NGF

储 存：冷冻（-20℃）

宿 主：Rabbit

相关类别：一抗

反应种属：Human, Mouse, Rat

标记物：Unconjugate

克隆类型：rabbit polyclonal

技术规格

Background:	This gene is a member of the NGF-beta family and encodes a secreted protein which homodimerizes and is incorporated into a larger complex. This protein has nerve growth stimulating activity and the complex is involved in the regulation of growth and the differentiation of sympathetic and certain sensory neurons. Mutations in this gene have been associated with hereditary sensory and autonomic neuropathy, type 5 (HSAN5), and dysregulation of this gene's expression is associated with allergic rhinitis.
Applications:	ELISA, IHC
Name of antibody:	NGF

Immunogen:	Fusion protein of human NGF
Full name:	nerve growth factor (beta polypeptide)
Synonyms :	NGFB; HSAN5; Beta-NGF
SwissProt:	P01138
ELISA Recommended dilution:	1000-5000
IHC positive control:	Human liver cancer
IHC Recommend dilution:	50-200

