

兔抗 NLRC4 多克隆抗体

中文名称：兔抗 NLRC4 多克隆抗体

英文名称：Anti-NLRC4 rabbit polyclonal antibody

别 名：NLR family CARD domain containing 4; CLAN; IPAF; AIFEC; CLAN1; CLANA; CLANB; CLANC; CLAND; FCAS4; CARD12; CLR2.1

相关类别：一抗

储 存：冷冻（-20℃）

宿 主：Rabbit

抗 原：NLRC4

反应种属：Human

标记物：Unconjugate

克隆类型：rabbit polyclonal

技术规格

Background:	This gene encodes a member of the caspase recruitment domain-containing NLR family. Family members play essential roles in innate immune response to a wide range of pathogenic organisms, tissue damage and other cellular stresses. Mutations in this gene result in autoinflammation with infantile enterocolitis. Alternative splicing results in multiple transcript variants.
Applications:	ELISA, IHC
Name of antibody:	NLRC4
Immunogen:	Fusion protein of human NLRC4

Full name:	NLR family CARD domain containing 4
Synonyms:	CLAN; IPAF; AIFEC; CLAN1; CLANA; CLANB; CLANC; CLAND; FCAS4; CARD12; CLR2.1
SwissProt:	Q9NPP4
ELISA Recommended dilution:	5000-10000
IHC positive control:	Human liver cancer
IHC Recommend dilution:	25-50

