

兔抗人 IgG

中文名称：兔抗人 IgG

英文名称：Rabbit anti-Human IgG

相关类别：其他二抗

储 存：冷冻 (-20℃) 避光

标 记 物：Unconjugate

技术规格

Background:	Immunoglobulin G (IgG) is one of the most abundant proteins in serum with normal levels between 8-17 mg/mL in adult blood. Ig G is important for our defense against microorganisms and the molecules are produced by B lymphocytes as a part of our adaptive immune response. The IgG molecule has two separate functions, to bind to the pathogen that elicited the response and to recruit other cells and molecules to destroy the antigen. The variability of the IgG pool is generated by somatic recombination and the number of specificities in an individual at a given time point is estimated to be 10 ¹¹ variants.
Applications:	1/1000-5000 for ELISA 1/1000-5000 for Western blotting 1/200-500 for Immunohistochemistry/cytochemistry