

兔抗 RFK 多克隆抗体

中文名称：兔抗 RFK 多克隆抗体

英文名称：Anti-RFK rabbit polyclonal antibody

别名：RIFK

相关类别：一抗

抗原：RFK

储存：冷冻（-20℃）

宿主：Rabbit

反应种属：Human, Mouse

标记物：Unconjugate

克隆类型：rabbit polyclonal

Background:	Riboflavin kinase (RFK; EC 2.7.1.26) is an essential enzyme that catalyzes the phosphorylation of riboflavin (vitamin B2) to form flavin mononucleotide (FMN), an obligatory step in vitamin B2 utilization and flavin cofactor synthesis (Karthikeyan et al., 2003 [PubMed 12623014]).
Applications:	ELISA, IHC
Name of antibody:	RFK
Immunogen:	Full length fusion protein
Full name:	riboflavin kinase
Synonyms :	RIFK
SwissProt:	Q969G6
ELISA Recommended dilution:	5000-10000

IHC positive control:	Human brain and Human liver cancer
IHC Recommend dilution:	25-100

