

兔抗 IGFBP4 多克隆抗体

中文名称：兔抗 IGFBP4 多克隆抗体

英文名称：Anti-IGFBP4 rabbit polyclonal antibody

别名：BP-4, IBP4, IGFBP-4, HT29-IGFBP

抗原：IGFBP4

储存：冷冻（-20℃）

宿主：Rabbit

相关类别：一抗

反应种属：Human, Mouse, Rat

标记物：Unconjugate

克隆类型：rabbit polyclonal

技术规格

IHC Recommend dilution:	25-100
IHC positive control:	Human liver cancer and human breast cancer
ELISA Recommended dilution:	2000-5000
SwissProt:	P22692
Synonyms :	BP-4; IBP4; IGFBP-4; HT29-IGFBP
Full name:	insulin-like growth factor binding protein 4
Immunogen:	Fusion protein of human IGFBP4
Name of antibody:	IGFBP4
Applications:	ELISA, IHC

Background:

This gene is a member of the insulin-like growth factor binding protein (IGFBP) family and encodes a protein with an IGFBP domain and a thyroglobulin type-I domain. The protein binds both insulin-like growth factors (IGFs) I and II and circulates in the plasma in both glycosylated and non-glycosylated forms. Binding of this protein prolongs the half-life of the IGFs and alters their interaction with cell surface receptors.



