

Cy3 标记的驴抗山羊 IgG

中文名称: Cy3 标记的驴抗山羊 IgG

英文名称: Cy3-conjugated Donkey Anti-Goat IgG

相关类别: 荧光标记二抗

储 存: 冷冻(-20℃) 避光

标 记 物: Cy3

技术规格

Applications	1/100-800 for most application
	rabbit and rat serum proteins, but it may cross-react with immunog lobulins from other species.
Background	tion with chicken, guinea pig, Syrian hamster, horse, human, mouse,
	by ELISA and/or solid-phase adsorbed to ensure minimal cross-reac
	on- immunoglobulin serum proteins. The antibody has been tested
	f other goat immunoglobulins. No antibody was detected against n
	with whole molecule goat IgG. It also reacts with the light chains o
	Based on immunoelectrophoresis and/or ELISA, the antibody reacts