

## Cy3 标记驴抗兔 IgG

中文名称: Cy3 标记驴抗兔 IgG

英文名称: Cy3-conjugated Donkey Anti-Rabbit IgG

相关类别: 荧光标记二抗

储 存: 冷冻 (-20℃)

标 记 物: Cy3

### 技术规格

<b>Background</b>	Based on immunoelectrophoresis and/or ELISA, the antibody reacts with whole molecule rabbit IgG. It also reacts with the light chains of other rabbit immunoglobulins. No antibody was detected against non-immunoglobulin serum proteins. The antibody has been tested by ELISA and/or solid-phase adsorbed to ensure minimal cross-reaction with bovine, chicken, goat, guinea pig, Syrian hamster, horse, human, mouse, rat, and sheep serum proteins, but it may cross-react with immunoglobulins from other species.
<b>Applications</b>	1/100-800 for most application