

## Cy3 标记驴抗小鼠 IgG

中文名称: Cy3 标记驴抗小鼠 IgG

英文名称: Cy3-conjugated Donkey Anti-Mouse IgG

相关类别: 荧光标记二抗

储 存: 冷冻(-20℃)

标 记 物: Cy3

## 技术规格

Background	Based on immunoelectrophoresis and/or ELISA, the antibody reacts with whole molecule mouse IgG. It also reacts with the light chains of other mouse immunoglobulins. No antibody was detected against non- immunoglobulin serum proteins. The antibody has been tested by ELISA and/or solid-phase adsorbed to ensure minimal cross-reacti on with bovine, chicken, goat, guinea pig, Syrian hamster, horse, hu man, rabbit, rat, and sheep serum proteins, but it may cross-react w ith immunoglobulins from other species.
Applications	1/100-800 for most application