

## Cy3 标记山羊抗鸡 IgG

- 中文名称: Cy3标记山羊抗鸡 IgG
- 英文名称: Cy3-conjugated Goat anti-chicken IgG
- 相关类别: 荧光标记二抗
- 储 存: 冷冻 (-20℃)
- 标记物: Cy3

## 技术规格

Background	Immunoglobulin G (IgG) is one of the most abundant proteins in se rum with normal levels between 8-17 mg/ml in adult blood. IgG is i mportant for our defense against microorganisms and the molecules are produced by B lymphocytes as a part of our adaptive immune r esponse. The IgG molecule has two separate functions; to bind to t he pathogen that elicited the response and to recruit other cells an d molecules to destroy the antigen. The variability of the IgG pool i s generated by somatic recombination and the number of specificiti es in an individual at a given time point is estimated to be 1011 va
	es in an individual at a given time point is estimated to be 1011 va riants.
Applications	1/100-1000 for FCM 1/100-1000 for Immunofluorescence