

兔抗 RGS19 多克隆抗体

中文名称：兔抗 RGS19 多克隆抗体

英文名称：Anti-RGS19 rabbit polyclonal antibody

别名：regulator of G-protein signaling 19 GAIP; RGSGAIP

相关类别：一抗

储存：冷冻（-20℃）

宿主：Rabbit

抗原：RGS19

反应种属：Human, Mouse, Rat

标记物：Unconjugate

克隆类型：rabbit polyclonal

技术规格

Background:	G proteins mediate a number of cellular processes. The protein encoded by this gene belongs to the RGS (regulators of G-protein signaling) family and specifically interacts with G protein, GAI3. This protein is a guanosine triphosphatase-activating protein that functions to down-regulate Galpha i/Galpha q-linked signaling. Alternatively spliced transcript variants encoding the same protein isoform have been found for this gene.
Applications:	ELISA, IHC
Name of antibody:	RGS19

Immunogen:	Synthetic peptide of human RGS19
Full name:	regulator of G-protein signaling 19
Synonyms:	GAIP; RGSGAIP
SwissProt:	P49795
ELISA Recommended dilution:	5000-10000
IHC positive control:	Human liver cancer;
IHC Recommend dilution:	25-100



