

兔抗 IGFBPL1 多克隆抗体

- 中文名称：兔抗 IGFBPL1 多克隆抗体
- 英文名称：Anti-IGFBPL1 rabbit polyclonal antibody
- 别名：insulin like growth factor binding protein like 1; IGFBPRP4; IGFBP-RP4; bA113O24.1
- 相关类别：一抗
- 储存：冷冻（-20℃）
- 宿主：Rabbit
- 抗原：IGFBPL1
- 反应种属：Human
- 标记物：Unconjugate
- 克隆类型：rabbit polyclonal

技术规格

Background:	IGF-binding proteins prolong the half-life of IGFs and have been shown to either inhibit or stimulate the growth promoting effects of the IGFs in cell culture. They alter the interaction of IGFs with their cell surface receptors (By similarity). May be a putative tumor suppressor protein.
Applications:	ELISA, IHC
Name of antibody:	IGFBPL1
Immunogen:	Synthetic peptide of human IGFBPL1
Full name:	insulin like growth factor binding protein like 1
Synonyms:	IGFBPRP4; IGFBP-RP4; bA113O24.1
SwissProt:	Q8WX77

ELISA Recommended dilution:	5000-10000
IHC positive control:	Human thyroid cancer
IHC Recommend dilution:	50-300

