

Cy3 标记的山羊抗豚鼠 IgG

- 中文名称: Cy3标记的山羊抗豚鼠 lgG
- 英文名称: Cy3-conjugated Goat anti-Guinea pig IgG
- 相关类别: 荧光标记二抗
- 储 存: 冷冻 (-20℃)
- 标记物: Cy3

技术规格

Background	Immunoglobulin G (IgG) is one of the most abundant proteins in serum with normal levels between 8-17 mg/ml in adult blood. IgG is important for our defense against microorganisms and the molecules are produced by B lymphocytes as a part of our adaptive immune response. The IgG molecule has two separate functions; to bind to the pathogen that elicit ed the response and to recruit other cells and molecules to destroy the antigen. The variability of the IgG pool is generated by somatic recombi nation and the number of specificities in an individual at a given time p oint is estimated to be 1011 variants.
Applications	1/20-100 for FCM 1/20-200 for Immunohistochemistry/cytochemistry