

Cy5.5 标记山羊抗人 IgG

- 中文名称: Cy5.5标记山羊抗人 IgG
- 英文名称: Cy5.5-conjugated Goat anti-human IgG
- 相关类别: 荧光标记二抗
- 储 存: 冷冻 (-20℃)
- 标记物: Cy5.5

技术规格

Background	Immunoglobulin G (IgG) is one of the most abundant proteins in serum with normal levels between 8-17 mg/ml in adult blood. Ig G is important for our defense against microorganisms and the m olecules are produced by B lymphocytes as a part of our adaptive immune response. The IgG molecule has two separate functions; t o bind to the pathogen that elicited the response and to recruit other cells and molecules to destroy the antigen. The variability o f the IgG pool is generated by somatic recombination and the nu mber of specificities in an individual at a given time point is esti mated to be 1011 variants.
Applications	1/20-100 for FCM 1/20-200 for Immunohistochemistry/cytochemist ry

