

Cy5 标记驴抗羊

中文名称: Cy5 标记驴抗羊

英文名称: Cy5-conjugated Donkey Anti-Goat IgG

相关类别: 荧光标记二抗

储 存: 冷冻(-20℃)

标 记 物: Cy5

技术规格

Background:	Immunoglobulin G (IgG) is one of the most abundant proteins in s erum with normal levels between 8-17 mg/mL in adult blood. IgG i s important for our defense against microorganisms and the molec ules are produced by B lymphocytes as a part of our adaptive imm une response. The IgG molecule has two separate functions, to bin d to the pathogen that elicited the response and to recruit other c ells and molecules to destroy the antigen. The variability of the Ig G pool is generated by somatic recombination and the number of specificities in an individual at a given time point is estimated to b e 1011 variants.
Applications:	1/100-1000 for Flow-Cyt and Immunofluorescence