

Cy5 标记山羊抗人 IgG

中文名称: Cy5 标记山羊抗人 IgG

英文名称: Cy5-conjugated Goat anti-human IgG

相关类别: 荧光标记二抗

储 存: 冷冻 (-20℃)

标 记 物: Cy5

技术规格

Background	Immunoglobulin G (IgG) is one of the most abundant proteins in serum with normal levels between 8-17 mg/ml in adult blood. Ig G is important for our defense against microorganisms and the molecules are produced by B lymphocytes as a part of our adaptive immune response. The IgG molecule has two separate functions; to bind to the pathogen that elicited the response and to recruit other cells and molecules to destroy the antigen. The variability of the IgG pool is generated by somatic recombination and the number of specificities in an individual at a given time point is estimated to be 10 ¹¹ variants.
Applications	1/20-100 for FCM 1/20-200 for Immunohistochemistry/cytochemistry

