

Cy5 标记山羊抗兔 IgG

中文名称: Cy5 标记山羊抗兔 IgG

英文名称: Cy5-conjugated Goat anti-Rabbit IgG

相关类别: 荧光标记二抗

储 存: 冷冻 (-20℃)

标 记 物: Cy5

技术规格

Background	Antibody Specificity: Based on immunoelectrophoresis and/or ELISA, the antibody reacts with whole molecule rabbit IgG. It also reacts with the light chains of other rabbit immunoglobulins. No antibody was detected against non-immunoglobulin serum proteins. The antibody has been tested by ELISA and/or solid-phase adsorption to ensure minimal cross-reaction with human, mouse, and rat serum proteins, but it may cross-react with immunoglobulins from other species.
Applications	1/50-200 for most application