

FITC 标记的驴抗绵羊 IgG

中文名称: FITC 标记的驴抗绵羊 IgG

英文名称: FITC-conjugated Donkey Anti-Sheep IgG

相关类别: 荧光标记二抗

储 存: 冷冻(-20℃)

标记物: FITC

技术规格

Background	Based on immunoelectrophoresis and/or ELISA, the antibody reacts with whole molecule sheep IgG. It also reacts with the light c hains of other sheep immunoglobulins. No antibody was detected against non-immunoglobulin serum proteins. The antibody has been tested by ELISA and/or solid-phase adsorbed to ensure minimal cross-reaction with chicken, guinea pig, Syrian hamster, horse, human, mouse, rabbit, rat serum proteins, but it may cross-react with immunoglobulins from other species.
Applications	1/50-200 for most application

