

FITC 标记的驴抗山羊 IgG

中文名称: FITC 标记的驴抗山羊 IgG

英文名称: FITC-conjugated Donkey Anti-Goat IgG

相关类别: 荧光标记二抗

储 存: 冷冻(-20℃) 避光

标记物: FITC

技术规格

Background	Based on immunoelectrophoresis and/or ELISA, the antibody r eacts with whole molecule goat IgG. It also reacts with the li ght chains of other goat immunoglobulins. No antibody was detected against non- immunoglobulin serum proteins. The a ntibody has been tested by ELISA and/or solid-phase adsorbe d to ensure minimal cross-reaction with chicken, guinea pig, Syrian hamster, horse, human, mouse, rabbit and rat serum p roteins, but it may cross-react with immunoglobulins from ot her species.
Applications	1/50-200 for most application