

兔抗 NPBWR1 多克隆抗体

中文名称：兔抗 NPBWR1 多克隆抗体

英文名称：Anti-NPBWR1 rabbit polyclonal antibody

别名：neuropeptides B and W receptor 1; GPR7

相关类别：一抗

储存：冷冻（-20℃）

宿主：Rabbit

抗原：NPBWR1

反应种属：Human, Mouse, Rat

标记物：Unconjugate

克隆类型：rabbit polyclonal

技术规格

Background:	Interacts specifically with a number of opioid ligands. Receptor for neuropeptides B and W, which may be involved in neuroendocrine system regulation, food intake and the organization of other signals. Has a higher affinity for neuropeptide B.
Applications:	ELISA, IHC
Name of antibody:	NPBWR1
Immunogen:	Synthetic peptide of human NPBWR1
Full name:	neuropeptides B and W receptor 1
Synonyms:	GPR7
SwissProt:	P48145

ELISA Recommended dilution:	5000-10000
IHC positive control:	Human tonsil
IHC Recommend dilution:	30-150

