

## FITC 标记山羊抗鸡 IgG

中文名称: FITC 标记山羊抗鸡 IgG

英文名称: FITC-conjugated Goat anti-chicken IgG

相关类别: 荧光标记二抗

储 存: 冷冻(-20℃)

标 记 物: FITC

## 技术规格

Background	Immunoglobulin G (IgG) is one of the most abundant proteins in se rum with normal levels between 8-17 mg/ml in adult blood. IgG is i mportant for our defense against microorganisms and the molecules are produced by B lymphocytes as a part of our adaptive immune r esponse. The IgG molecule has two separate functions; to bind to the pathogen that elicited the response and to recruit other cells and molecules to destroy the antigen. The variability of the IgG pool is generated by somatic recombination and the number of specificities in an individual at a given time point is estimated to be 1011 variants.
Applications	1/100-1000 for Flow-Cyt and Immunofluorescence