

## FITC 标记山羊抗人 IgG

中文名称： FITC 标记山羊抗人 IgG

英文名称： FITC-conjugated Goat anti-human IgG

相关类别： 荧光标记二抗

储 存 冷冻（-20℃）

标 记 物： FITC

### 技术规格

<b>Background</b>	Antibody Specificity: Based on immunoelectrophoresis and/or ELISA, the antibody reacts with whole molecule human IgG. It also reacts with the light chains of other human immunoglobulins. No antibody was detected against non-immunoglobulin serum proteins. The antibody may cross-react with immunoglobulins from other species.
<b>Applications</b>	1/50-200 for most application