

兔抗 IKBKE 多克隆抗体

- 中文名称：兔抗 IKBKE 多克隆抗体
- 英文名称：Anti-ICKBE rabbit polyclonal antibody
- 别名：IKKE; IKKI; IKK-E; IKK-i
- 相关类别：一抗
- 储存：冷冻（-20℃）
- 抗原：IKBKE
- 宿主：Rabbit
- 反应种属：Human, Mouse
- 标记物：Unconjugate
- 克隆类型：rabbit polyclonal

技术规格

Background:	IKBKE is a noncanonical I-kappa-B (see MIM 164008) kinase (I KK) that is essential for regulating antiviral signaling pathways. IKBKE has also been identified as a breast cancer (MIM 11448 0) oncogene and is amplified and overexpressed in over 30% of breast carcinomas and breast cancer cell lines
Applications:	ELISA, IHC
Name of antibody:	IKBKE
Immunogen:	Fusion protein of human IKBKE
Full name:	inhibitor of kappa light polypeptide gene enhancer in B-cells, kinase epsilon
Synonyms :	IKKE; IKKI; IKK-E; IKK-i

SwissProt:	Q14164
ELISA Recommended dilution:	2000-5000
IHC positive control:	Human ovarian cancer and Human thyroid cancer
IHC Recommend dilution:	50-200

