

FITC 标记兔抗狗 IgG

- 中文名称: FITC 标记兔抗狗 lgG
- 英文名称: FITC-conjugated Rabbit Anti-Dog IgG
- 相关类别: 荧光标记二抗
- 储 存: 冷冻 (-20℃)
- 标记物: FITC

技术规格

Background:	Immunoglobulin G (IgG) is one of the most abundant proteins i n serum with normal levels between 8-17 mg/mL in adult blood . IgG is important for our defense against microorganisms and t he molecules are produced by B lymphocytes as a part of our a daptive immune response. The IgG molecule has two separate f unctions, to bind to the pathogen that elicited the response an d to recruit other cells and molecules to destroy the antigen. T he variability of the IgG pool is generated by somatic recombin ation and the number of specificities in an individual at a given time point is estimated to be 1011 variants.
Applications:	1/100-1000 for FCM 1/100-1000 for Immunofluorescence