

FITC 标记小鼠抗山羊 IgG

中文名称: FITC 标记小鼠抗山羊 IgG

英文名称: FITC-conjugated Mouse anti-goat IgG

相关类别: 荧光标记二抗

储 存: 冷冻 (-20℃)

标 记 物: FITC

技术规格

Background	Antibody Specificity: Based on immunoelectrophoresis and/or ELISA, the antibody reacts with whole molecule goat IgG. It also reacts with the light chains of other goat immunoglobulins. No antibody was detected against non-immunoglobulin serum proteins. The antibody exhibits inherent minimal cross-reaction to mouse serum proteins and has been tested by ELISA and/or solid-phase adsorbed to ensure minimal cross-reaction with mouse, human, and rabbit serum proteins. The antibody may cross-react with immunoglobulins from other species.
Applications	1/50-200 for most application