

RBITC 标记山羊抗兔 IgG

中文名称: RBITC 标记山羊抗兔 IgG

英文名称: RBITC-conjugated Goat anti-Rabbit IgG

相关类别: 荧光标记二抗

储 存: 冷冻(-20℃)

标记物: RBITC

技术规格

Background	Immunoglobulin G (IgG) is one of the most abundant proteins in ser um with normal levels between 8-17 mg/ml in adult blood. IgG is i mportant for our defense against microorganisms and the molecules are produced by B lymphocytes as a part of our adaptive immune r esponse. The IgG molecule has two separate functions; to bind to the pathogen that elicited the response and to recruit other cells and molecules to destroy the antigen. The variability of the IgG pool is g enerated by somatic recombination and the number of specificities i n an individual at a given time point is estimated to be 1011 variant s.
Applications	1/500-1000 for ELISA 1/20-200 for Immunohistochemistry/Immunocyt ocehmistry