

TRITC 标记驴抗鸡 IgG

中文名称: TRITC 标记驴抗鸡 IgG

英文名称: TRITC-conjugated Donkey Anti-Chicken IgG

相关类别: 荧光标记二抗

储 存: 冷冻(-20℃)

标 记 物: TRITC

技术规格

Background	Based on immunoelectrophoresis and/or ELISA, the antibody re acts with whole molecule Chicken IgY. It also reacts with the li ght chains of other Chicken immunoglobulins. No antibody was detected against non-immunoglobulin serum proteins. The antibody has been tested by ELISA and/or solid-phase adsorbed to ensure minimal cross-reaction with bovine, goat, guinea pig, Syrian hamster, horse, human, mouse, rabbit, rat, and sheep se rum proteins, but it may cross-react with immunoglobulins from other species.
Applications	1/50-200 for most application

