

兔抗 NQO2 多克隆抗体

中文名称：兔抗 NQO2 多克隆抗体

英文名称：Anti-NQO2 rabbit polyclonal antibody

别名：QR2, DHQV, DIA6, NMOR2

相关类别：一抗

储存：冷冻（-20℃）

宿主：Rabbit

抗原：NQO2

反应种属：Human, Mouse


标记物：Unconjugate



克隆类型：rabbit polyclonal

技术规格

Background:	NQO2 (EC 1.10.99.2) is a flavoprotein that catalyzes the 2-electron reduction of various quinones, redox dyes, and the vitamin K menadione. NQO2 predominantly uses dihydronicotinamide riboside (NRH) as the electron donor. The enzyme apparently serves as a quinone reductase in connection with conjugation reactions of hydroquinones involved in detoxification pathways as well as in biosynthetic processes such as the vitamin K-dependent gamma-carboxylation of glutamate residues in prothrombin synthesis.
Applications:	ELISA, WB, IHC
Name of antibody:	NQO2
Immunogen:	Synthetic peptide of human NQO2

Full name:	NAD(P)H dehydrogenase, quinone 2
Synonyms :	QR2, DHQV, DIA6, NMOR2
SwissProt:	P16083
ELISA Recommended dilution:	2000-10000
IHC positive control:	Human breast cancer and Human liver cancer
IHC Recommend dilution:	50-200
WB Predicted band size:	26 kDa
WB Positive control:	Mouse liver tissue
WB Recommended dilution:	1000-5000





mlbio 梅联生物
Good elisakit producers