

## 兔抗 NR2F2 多克隆抗体

中文名称：兔抗 NR2F2 多克隆抗体

英文名称： Anti-NR2F2 rabbit polyclonal antibody

别名： ARP1; CHTD4; NF-E3; NR2F1; SVP40; COUPTFB; TFCOUP2; COUPTFII

相关类别： 一抗

储存： 冷冻（-20℃）

宿主： Rabbit

抗原： NR2F2

反应种属： Human, Mouse, Rat

标记物： Unconjugate

克隆类型： rabbit polyclonal

技术规格

<b>Background:</b>	This gene encodes a member of the steroid thyroid hormone superfamily of nuclear receptors. The encoded protein is a ligand inducible transcription factor that is involved in the regulation of many different genes. Alternate splicing results in multiple transcript variants.
<b>Applications:</b>	ELISA, WB
<b>Name of antibody:</b>	NR2F2
<b>Immunogen:</b>	Synthetic peptide of humanNR2F2
<b>Full name:</b>	nuclear receptor subfamily 2, group F, member 2
<b>Synonyms :</b>	ARP1; CHTD4; NF-E3; NR2F1; SVP40; COUPTFB; TFCOUP2 ; COUPTFII

<b>SwissProt:</b>	P24468
<b>ELISA Recommended dilution:</b>	1000-2000
<b>WB Predicted band size:</b>	46 kDa
<b>WB Positive control:</b>	Human placenta tissue
<b>WB Recommended dilution:</b>	200-1000

