

兔抗 IL18 多克隆抗体

中文名称：兔抗 IL18 多克隆抗体

英文名称：Anti-IL18 rabbit polyclonal antibody

别名：IGIF, IL-18, IL-1g, IL1F4

抗原：IL18

储存：冷冻（-20℃）

宿主：Rabbit

相关类别：一抗

反应种属：Human

标记物：Unconjugate

克隆类型：rabbit polyclonal

技术规格

Background:	The protein encoded by this gene is a proinflammatory cytokine that augments natural killer cell activity in spleen cells, and stimulates interferon gamma production in T-helper type I cells. Alternatively spliced transcript variants encoding different isoforms have been found for this gene
Applications:	ELISA, WB, IHC
Name of antibody:	IL18
Immunogen:	Fusion protein of human IL18
Full name:	interleukin 18
Synonyms :	IGIF, IL-18, IL-1g, IL1F4
SwissProt:	Q14116

ELISA Recommended dilution:	2000-5000
IHC positive control:	Human ovarian cancer and Human cervical cancer
IHC Recommend dilution :	100-200
WB Predicted band size:	22 kDa
WB Positive control:	Hela cells
WB Recommended dilution:	1000-2000



