

## 兔抗 RPL12 多克隆抗体

中文名称：兔抗 RPL12 多克隆抗体

英文名称：Anti-RPL12 rabbit polyclonal antibody

别名：L12

抗原：RPL12

储存：冷冻（-20℃）避光

宿主：Rabbit

反应种属：Human Mouse

相关类别：一抗

标记物：Unconjugate

克隆类型：Unconjugate

### 技术规格

#### Background:

Ribosomes, the organelles that catalyze protein synthesis, consist of a small 40S subunit and a large 60S subunit. Together these subunits are composed of 4 RNA species and approximately 80 structurally distinct proteins. This gene encodes a ribosomal protein that is a component of the 60S subunit. The protein belongs to the L11P family of ribosomal proteins. It is located in the cytoplasm. The protein binds directly to the 26S rRNA. This gene is co-transcribed with the U65 snoRNA, which is located in its fourth intron. As is typical for genes encoding ribosomal protein

	s, there are multiple processed pseudogenes of this gene dispersed through the genome.
<b>Applications:</b>	WB
<b>Name of antibody:</b>	RPL12
<b>Immunogen:</b>	Synthesized peptide derived from internal of human RPL12.
<b>Full name:</b>	ribosomal protein L12
<b>Synonyms :</b>	L12
<b>SwissProt:</b>	P30050
<b>WB Predicted band size:</b>	18 kDa
<b>WB Positive control:</b>	HepG2 cells and Hela cells lysates
<b>WB Recommended dilution:</b>	500-3000

