

## 兔抗 RPL18 多克隆抗体

中文名称：兔抗 RPL18 多克隆抗体

英文名称：Anti-RPL18 rabbit polyclonal antibody

别名：ribosomal protein L18; L18

抗原：RPL18

储存：冷冻（-20℃）

宿主：Rabbit

反应种属：Human, Mouse, Rat

相关类别：一抗

标记物：Unconjugate

克隆类型：rabbit polyclonal

### 技术规格

#### Background:

Ribosomes, the organelles that catalyze protein synthesis, consist of a small 40S subunit and a large 60S subunit. Together these subunits are composed of 4 RNA species and approximately 80 structurally distinct proteins. This gene encodes a member of the L18E family of ribosomal proteins that is a component of the 60S subunit. As is typical for genes encoding ribosomal proteins, there are multiple processed pseudogenes of this gene dispersed through the genome. Alternatively spliced transcript variants encoding multiple isoforms have been observed for this gene.

<b>Applications:</b>	ELISA, IHC
<b>Name of antibody:</b>	RPL18
<b>Immunogen:</b>	Fusion protein of human RPL18
<b>Full name:</b>	ribosomal protein L18
<b>Synonyms:</b>	L18
<b>SwissProt:</b>	Q07020
<b>ELISA Recommended dilution:</b>	5000-10000
<b>IHC positive control:</b>	Human liver cancer and Human colorectal cancer
<b>IHC Recommend dilution:</b>	50-200



