

## 兔抗 NTRK2 多克隆抗体

中文名称：兔抗 NTRK2 多克隆抗体

英文名称：Anti-NTRK2 rabbit polyclonal antibody

别名：neurotrophic receptor tyrosine kinase 2 TRKB; trk-B; GP145-TrkB

相关类别：一抗

储存：冷冻（-20℃）

宿主：Rabbit

抗原：NTRK2

反应种属：Human, Mouse, Rat

标记物：Unconjugate

克隆类型：rabbit polyclonal

### 技术规格

<b>Background:</b>	This gene encodes a member of the neurotrophic tyrosine receptor kinase (NTRK) family. This kinase is a membrane-bound receptor that, upon neurotrophin binding, phosphorylates itself and members of the MAPK pathway. Signaling through this kinase leads to cell differentiation. Mutations in this gene have been associated with obesity and mood disorders. Alternative splicing results in multiple transcript variants.
<b>Applications:</b>	ELISA, IHC
<b>Name of antibody:</b>	NTRK2
<b>Immunogen:</b>	Fusion protein of human NTRK2
<b>Full name:</b>	neurotrophic receptor tyrosine kinase 2
<b>Synonyms:</b>	TRKB; trk-B; GP145-TrkB

<b>SwissProt:</b>	Q16620
<b>ELISA Recommended dilution:</b>	5000-10000
<b>IHC positive control:</b>	Human liver cancer
<b>IHC Recommend dilution:</b>	20-100

