

兔抗 UQCRQ 多克隆抗体

中文名称：兔抗 UQCRQ 多克隆抗体

英文名称：Anti-UQCRQ rabbit polyclonal antibody

别名：QPC; QCR8; QP-C; UQCR7; MC3DN4

相关类别：一抗

抗原：UQCRQ

储存：冷冻（-20℃）

宿主：Rabbit

反应种属：Human, Rat

标记物：Unconjugate

克隆类型：rabbit polyclonal

技术规格

Background:	This gene encodes a ubiquinone-binding protein of low molecular mass. This protein is a small core-associated protein and a subunit of ubiquinol-cytochrome c reductase complex III, which is part of the mitochondrial respiratory chain.
Applications:	ELISA, WB, IHC
Name of antibody:	UQCRQ
Immunogen:	Full length fusion protein
Full name:	ubiquinol-cytochrome c reductase complex III subunit VII
Synonyms:	QPC; QCR8; QP-C; UQCR7; MC3DN4
SwissProt:	O14949

ELISA Recommended dilution:	5000-10000
IHC positive control:	Human liver cancer
IHC Recommend dilution:	40-200
WB Predicted band size:	10 kDa
WB Positive control:	293T , A549 , Hela and 231 cell , Human placenta tissue , HEPG2 cell , Human fetal liver tissue lysates
WB Recommended dilution:	500-2000



