

兔抗 RPS14 多克隆抗体

中文名称: 兔抗 RPS14 多克隆抗体

英文名称: Anti-RPS14 rabbit polyclonal antibody

别 名: ribosomal protein S14; S14; EMTB

相关类别: 一抗

储 存: 冷冻(-20℃)

宿 主: Rabbit

抗 原: RPS14

反应种属: Human, Mouse, Rat

标记物: Unconjugate

克隆类型: rabbit polyclonal

技术规格

Background:

Ribosomes, the organelles that catalyze protein synthesis, consist of a small 40S subunit and a large 60S subunit. Together these subunits are composed of 4 RNA species and approximately 80 structurally distinct proteins. This gene encodes a ribosomal protein that is a component of the 40S subunit. The protein belongs to the S1 1P family of ribosomal proteins. It is located in the cytoplasm. Transcript variants utilizing alternative transcription initiation sites have been described in the literature. As is typical for genes encoding ribosomal proteins, the re are multiple processed pseudogenes of this gene dis



	persed through the genome. In Chinese hamster ovary cells, mutations in this gene can lead to resistance to e metine, a protein synthesis inhibitor. Multiple alternative ly spliced transcript variants encoding the same protein have been found for this gene.
Applications:	ELISA, WB, IHC
Name of antibody:	RPS14
Immunogen:	Fusion protein of human RPS14
Full name:	ribosomal protein S14
Synonyms:	S14; EMTB
SwissProt:	P62263
ELISA Recommended dilution:	5000-10000
IHC positive control:	Human liver cancer and Human cervical cancer
IHC Recommend dilution:	50-200
WB Predicted band size:	16 kDa
WB Positive control:	231, HL60 and Jurkat cell lysates
WB Recommended dilution:	1000-5000







