

兔抗 ODC1 多克隆抗体

中文名称：兔抗 ODC1 多克隆抗体

英文名称：Anti-ODC1 rabbit polyclonal antibody

别名：ODC

抗原：ODC1

储存：冷冻（-20℃）

宿主：Rabbit

相关类别：一抗

反应种属：Human, Mouse, Rat

标记物：Unconjugate

克隆类型：rabbit polyclonal

技术规格

Background:	This gene encodes the rate-limiting enzyme of the polyamine biosynthesis pathway which catalyzes ornithine to putrescine. The activity level for the enzyme varies in response to growth-promoting stimuli and exhibits a high turnover rate in comparison to other mammalian proteins. Originally localized to both chromosomes 2 and 7, the gene encoding this enzyme has been determined to be located on 2p25, with a pseudogene located on 7q31-qter.
Applications:	ELISA, IHC
Name of antibody:	ODC1
Immunogen:	Fusion protein of human ODC1
Full name:	ornithine decarboxylase 1
Synonyms :	ODC

SwissProt:	P11926
ELISA Recommended dilution:	1000-5000
IHC positive control:	Human colon cancer and Human thyroid cancer
IHC Recommend dilution:	25-100

