

## 兔抗 USP50 多克隆抗体

中文名称：兔抗 USP50 多克隆抗体

英文名称：Anti-USP50 rabbit polyclonal antibody

别名：ubiquitin specific peptidase 50

相关类别：一抗

储存：冷冻（-20℃）

宿主：Rabbit

抗原：USP50

反应种属：Human, Mouse

标记物：Unconjugate

克隆类型：rabbit polyclonal

### 技术规格

<b>Background:</b>	Has no peptidase activity.
<b>Applications:</b>	ELISA, WB
<b>Name of antibody:</b>	USP50
<b>Immunogen:</b>	Synthetic peptide of human USP50
<b>Full name:</b>	ubiquitin specific peptidase 50
<b>SwissProt:</b>	Q70EL3
<b>ELISA Recommended dilution:</b>	5000-10000
<b>WB Predicted band size:</b>	39 kDa
<b>WB Positive control:</b>	Mouse adrenal gland tissue lysate
<b>WB Recommended dilution:</b>	500-2000

