

兔抗 INTS11 多克隆抗体

- 中文名称： 兔抗 INTS11 多克隆抗体
- 英文名称： Anti-INTS11 rabbit polyclonal antibody
- 别 名： integrator complex subunit 11; RC68; INT11; RC-68;
- 相关类别： 一抗
- 储 存： 冷冻（-20℃）
- 宿 主： Rabbit
- 抗 原： INTS11
- 反应种属： Human, Mouse, Rat
- 标记物： Unconjugate
- 克隆类型： rabbit polyclonal

技术规格

Background:	The Integrator complex contains at least 12 subunits and associates with the C-terminal domain of RNA polymerase II large subunit (POLR2A; MIM 180660) and mediates the 3-prime end processing of small nuclear RNAs U1 (RNU1; MIM 180680) and U2 (RNU2; MIM 180690). INTS11, or CPSF3L, is the catalytic subunit of the Integrator complex (Baillat et al., 2005 [PubMed 16239144]).[supplied by OMIM, Mar 2008]
Applications:	ELISA, IHC
Name of antibody:	INTS11

Immunogen:	Fusion protein of human INTS11
Full name:	integrator complex subunit 11
Synonyms:	RC68; INT11; RC-68; CPSF3L; CPSF73L
SwissProt:	Q5TA45
IHC positive control:	Human ovarian cancer
IHC Recommend dilution:	25-100

