

## 兔抗 UTP3 多克隆抗体

中文名称：兔抗 UTP3 多克隆抗体

英文名称：Anti-UTP3 rabbit polyclonal antibody

别名：CRL1; CRLZ1; SAS10

相关类别：一抗

储存：冷冻（-20℃）

抗原：UTP3

宿主：Rabbit

反应种属：Human, Mouse, Rat

标记物：Unconjugate

克隆类型：rabbit polyclonal

### 技术规格

|                                    |                                                                                                                                        |
|------------------------------------|----------------------------------------------------------------------------------------------------------------------------------------|
| <b>Background:</b>                 | Essential for gene silencing: has a role in the structure of silenced chromatin. Plays a role in the developing brain (By similarity). |
| <b>Applications:</b>               | ELISA, IHC                                                                                                                             |
| <b>Name of antibody:</b>           | UTP3                                                                                                                                   |
| <b>Immunogen:</b>                  | Fusion protein of human UTP3                                                                                                           |
| <b>Full name:</b>                  | UTP3 small subunit processome component                                                                                                |
| <b>Synonyms:</b>                   | CRL1; CRLZ1; SAS10                                                                                                                     |
| <b>SwissProt:</b>                  | Q9NQZ2                                                                                                                                 |
| <b>ELISA Recommended dilution:</b> | 5000-10000                                                                                                                             |
| <b>IHC positive control:</b>       | Human colorectal cancer and Human tonsil                                                                                               |

IHC Recommend dilution:

40-200



