

## 兔抗 ORM1 多克隆抗体

中文名称：兔抗 ORM1 多克隆抗体

英文名称：Anti-ORM1 rabbit polyclonal antibody

别名：ORM; AGP1; AGP-A

相关类别：一抗

储存：冷冻（-20℃）

抗原：ORM1

宿主：Rabbit

反应种属：Human

标记物：Unconjugate

克隆类型：rabbit polyclonal

### 技术规格

|                          |   |
|--------------------------|---|
| <b>Background:</b>       | This gene encodes a key acute phase plasma protein. Because of its increase due to acute inflammation, this protein is classified as an acute-phase reactant. The specific function of this protein has not yet been determined; however, it may be involved in aspects of immunosuppression. |
| <b>Applications:</b>     | ELISA, WB   |
| <b>Name of antibody:</b> | ORM1  |
| <b>Immunogen:</b>        | Fusion protein of human ORM1  |
| <b>Full name:</b>        | orosomucoid 1   |
| <b>Synonyms :</b>        | ORM; AGP1; AGP-A  |
| <b>SwissProt:</b>        | P02763  |

|                                    |                                  |
|------------------------------------|----------------------------------|
| <b>ELISA Recommended dilution:</b> | 1000-2000                        |
| <b>WB Predicted band size:</b>     | 24 kDa                           |
| <b>WB Positive control:</b>        | Human testis and prostate tissue |
| <b>WB Recommended dilution:</b>    | 200-1000                         |

