

兔抗 IRAK1 多克隆抗体

中文名称：兔抗 IRAK1 多克隆抗体

英文名称：Anti-IRAK1 rabbit polyclonal antibody

别名：IRAK; pelle

相关类别：一抗

储存：冷冻（-20℃）

宿主：Rabbit

抗原：IRAK1

反应种属：Human

标记物：Unconjugate

克隆类型：rabbit polyclonal

技术规格

Background:	This gene encodes the interleukin-1 receptor-associated kinase 1, one of two putative serine/threonine kinases that become associated with the interleukin-1 receptor (IL1R) upon stimulation. This gene is partially responsible for IL1-induced upregulation of the transcription factor NF-kappa B. Alternatively spliced transcript variants encoding different isoforms have been found for this gene.
Applications:	ELISA, IHC
Name of antibody:	IRAK1
Immunogen:	Synthetic peptide of human IRAK1
Full name:	interleukin 1 receptor associated kinase 1

Synonyms:	IRAK; pelle
SwissProt:	P51617
ELISA Recommended dilution:	5000-10000
IHC positive control:	Human liver cancer and human thyroid cancer
IHC Recommend dilution:	25-100



