

兔抗 OXSM 多克隆抗体

中文名称：兔抗 OXSM 多克隆抗体

英文名称：Anti-OXSM rabbit polyclonal antibody

别名：3-oxoacyl-ACP synthase, mitochondrial; KS; CEM1; KASI; FASN2D

相关类别：一抗

储存：冷冻（-20℃）

宿主：Rabbit

抗原：OXSM

反应种属：Human, Mouse

标记物：Unconjugate

克隆类型：rabbit polyclonal

技术规格

Background:	This gene encodes a beta-ketoacyl synthetase. The encoded enzyme is required for elongation of fatty acid chains in the mitochondria. Alternatively spliced transcript variants have been described.
Applications:	ELISA, WB, IHC
Name of antibody:	OXSM
Immunogen:	Fusion protein of human OXSM
Full name:	3-oxoacyl-ACP synthase, mitochondrial
Synonyms:	KS; CEM1; KASI; FASN2D
SwissProt:	Q9NWU1
ELISA Recommended dilution	5000-10000

:	
IHC positive control:	Human colorectal cancer and Human esophagus cancer
IHC Recommend dilution:	100-300
WB Predicted band size:	49 kDa
WB Positive control:	NIH/3T3 cell lysate
WB Recommended dilution:	200-1000



