

## 兔抗 OXSR1(Phospho-Thr185) 多克隆抗体

中文名称：兔抗 OXSR1(Phospho-Thr185) 多克隆抗体

英文名称：Anti-OXSR1(Phospho-Thr185) rabbit polyclonal antibody

别 名：OSR1

相关类别：一抗

储 存：冷冻（-20℃）避光

宿 主：Rabbit

抗 原：OXSR1(Phospho-Thr185)

反应种属：Human

标记物：Unconjugate

克隆类型：rabbit polyclonal

### 技术规格

<b>Background:</b>	The product of this gene belongs to the Ser/Thr protein kinase family of proteins. It regulates downstream kinases in response to environmental stress, and may play a role in regulating the actin cytoskeleton.
<b>Applications:</b>	WB, IHC
<b>Name of antibody:</b>	OXSR1(Phospho-Thr185)
<b>Immunogen:</b>	Peptide sequence around phosphorylation site of threonine 185(R-K-T(p)-F-V) derived from Human OSR1.
<b>Full name:</b>	oxidative stress responsive 1
<b>Synonyms :</b>	OSR1
<b>SwissProt:</b>	O95747

<b>IHC positive control:</b>	Human breast carcinoma tissue
<b>IHC Recommend dilution:</b>	50-100
<b>WB Predicted band size:</b>	58 kDa
<b>WB Positive control:</b>	HepG2 cells lysates treated with serum
<b>WB Recommended dilution:</b>	500-1000



