

兔抗 P3H3 多克隆抗体

中文名称：兔抗 P3H3 多克隆抗体

英文名称：Anti-P3H3 rabbit polyclonal antibody

别名：prolyl 3-hydroxylase 3; GRCB; LEPREL2; HSU47926

相关类别：一抗

储存：冷冻（-20℃）

抗原：P3H3

宿主：Rabbit

反应种属：Human, Mouse

标记物：Unconjugate

克隆类型：rabbit polyclonal

技术规格

Background:	The protein encoded by this gene belongs to the leprecan family of proteoglycans, which function as collagen prolyl hydroxylases that are required for proper collagen biosynthesis, folding and assembly. This protein, like other family members, is thought to reside in the endoplasmic reticulum. Epigenetic inactivation of this gene is associated with breast and other cancers, suggesting that it may function as a tumor suppressor. [provided by RefSeq, Aug 2013]
Applications:	ELISA, IHC
Name of antibody:	P3H3
Immunogen:	Fusion protein of human P3H3

Full name:	prolyl 3-hydroxylase 3
Synonyms:	GRCB; LEPREL2; HSU47926
SwissProt:	Q8IVL6
IHC positive control:	Human colon cancer and Human brain
IHC Recommend dilution:	50-200

