

兔抗 PAF1 多克隆抗体

中文名称：兔抗 PAF1 多克隆抗体

英文名称：Anti-PAF1 rabbit polyclonal antibody

别名：PAF1 homolog, Paf1/RNA polymerase II complex component; PD2; F23149_1

相关类别：一抗

储存：冷冻（-20℃）

宿主：Rabbit

抗原：PAF1

反应种属：Human, Mouse, Rat

标记物：Unconjugate

克隆类型：rabbit polyclonal

技术规格

Background:	This gene encodes a subunit of the polymerase associated factor (PAF1) complex. The PAF1 complex interacts with RNA polymerase II and plays a role in transcription elongation as well as histone modifications including ubiquitylation and methylation. Alternatively spliced transcript variants encoding multiple isoforms have been observed for this gene.
Applications:	ELISA, IHC
Name of antibody:	PAF1
Immunogen:	Fusion protein of human PAF1
Full name:	PAF1 homolog, Paf1/RNA polymerase II complex comp

	onent
Synonyms:	PD2; F23149_1
SwissProt:	Q8N7H5
ELISA Recommended dilution:	5000-10000
IHC positive control:	Human cervical cancer and Human thyroid cancer
IHC Recommend dilution:	50-300



