

兔抗 IRAK3 多克隆抗体

- 中文名称：兔抗 IRAK3 多克隆抗体
- 英文名称：Anti-IRAK3 rabbit polyclonal antibody
- 别名：ASRT5; IRAKM
- 抗原：IRAK3
- 储存：冷冻（-20℃）
- 宿主：Rabbit
- 相关类别：一抗
- 反应种属：Human, Mouse
- 标记物：Unconjugate
- 克隆类型：rabbit polyclonal

技术规格

Background:	This gene encodes a member of the interleukin-1 receptor-associated kinase protein family. Members of this family are essential components of the Toll/IL-R immune signal transduction pathways. This protein is primarily expressed in monocytes and macrophages and functions as a negative regulator of Toll-like receptor signaling. Mutations in this gene are associated with a susceptibility to asthma. Alternate splicing results in multiple transcript variants.
Applications:	ELISA, IHC
Name of antibody:	IRAK3
Immunogen:	Fusion protein of human IRAK3

Full name:	interleukin-1 receptor-associated kinase 3
Synonyms :	ASRT5; IRAKM
SwissProt:	Q9Y616
ELISA Recommended dilution:	1000-2000
IHC positive control:	Human thyroid cancer
IHC Recommend dilution:	25-100

