

兔抗 PAK2(Ab-141) 多克隆抗体

中文名称：兔抗 PAK2(Ab-141) 多克隆抗体

英文名称：Anti-PAK2(Ab-141) rabbit polyclonal antibody

别名：PAK65; PAKgamma

相关类别：一抗

储存：冷冻（-20℃）避光

抗原：PAK2(Ab-141)

宿主：Rabbit

反应种属：Human

标记物：Unconjugate

克隆类型：rabbit polyclonal

技术规格

Background:	The p21 activated kinases (PAK) are critical effectors that link Rho GTPases to cytoskeleton reorganization and nuclear signaling. The PAK proteins are a family of serine/threonine kinases that serve as targets for the small GTP binding proteins, CDC42 and RAC1, and have been implicated in a wide range of biological activities. The protein encoded by this gene is activated by proteolytic cleavage during caspase-mediated apoptosis, and may play a role in regulating the apoptotic events in the dying cell.
Applications:	WB
Name of antibody:	PAK2(Ab-141)

Immunogen:	Synthesized non-phosphopeptide derived from human PA K2 around the phosphorylation site of serine 141 (Y-L-S(p)-F-T).
Full name:	p21 protein (Cdc42/Rac)-activated kinase 2
Synonyms :	PAK65; PAKgamma
SwissProt:	Q13177
WB Predicted band size:	58 kDa
WB Positive control:	Hela cells treated with TSA (400nM, 24hours) lysate
WB Recommended dilution:	500-3000

