

兔抗 S100Z 多克隆抗体

中文名称: 兔抗 S100Z 多克隆抗体

英文名称: Anti-S100Z rabbit polyclonal antibody

别名: Gm625; S100-zeta

储存: 冷冻 (-20℃) 避光

抗原: S100Z

宿主: Rabbit

反应种属: Human

相关类别: 一抗

标记物: Unconjugate

克隆类型: Unconjugate

技术规格

Background:	Members of the S100 protein family contain 2 calcium-binding EF-hands and exhibit cell-type specific expression patterns. For additional background information on S100 proteins, see MIM 114085.
Applications:	WB
Name of antibody:	S100Z
Immunogen:	Synthesized peptide derived from internal of human S100Z.
Full name:	S100 calcium binding protein Z
Synonyms :	Gm625; S100-zeta
SwissProt:	Q8WVG8

WB Predicted band size:	12 kDa
WB Positive control:	COLO cells lysate
WB Recommended dilution:	500-3000

