

## 兔抗 SCAND1 多克隆抗体

中文名称：兔抗 SCAND1 多克隆抗体

英文名称：Anti-SCAND1 rabbit polyclonal antibody

别名：SCAN domain containing 1; RAZ1; SDP1

相关类别：一抗

储存：冷冻（-20℃）

宿主：Rabbit

抗原：SCAND1

反应种属：Human

标记物：Unconjugate

克隆类型：rabbit polyclonal

### 技术规格

<b>Background:</b>	This gene encodes a SCAN box domain-containing protein. The SCAN domain is a highly conserved, leucine-rich motif of approximately 60 aa originally found within a subfamily of zinc finger proteins. This gene belongs to a family of genes that encode an isolated SCAN domain, but no zinc finger motif. This protein binds to and may regulate the function of the transcription factor myeloid zinc finger 1B. Alternate splicing results in multiple transcript variants.
<b>Applications:</b>	ELISA, IHC
<b>Name of antibody:</b>	SCAND1

<b>Immunogen:</b>	Fusion protein of human SCAND1
<b>Full name:</b>	SCAN domain containing 1
<b>Synonyms:</b>	RAZ1; SDP1
<b>SwissProt:</b>	P57086
<b>ELISA Recommended dilution:</b>	5000-10000
<b>IHC positive control:</b>	Human esophagus cancer
<b>IHC Recommend dilution:</b>	50-300

