

兔抗 PDAP1 多克隆抗体

中文名称：兔抗 PDAP1 多克隆抗体

英文名称：Anti-PDAP1 rabbit polyclonal antibody

别名：PDGFA associated protein 1; PAP; PAP1; HASPP28

相关类别：一抗

储存：冷冻（-20℃）

宿主：Rabbit

抗原：PDAP1

反应种属：Human, Mouse, Rat

标记物：Unconjugate

克隆类型：rabbit polyclonal

技术规格

Background:	The protein encoded by this gene is a phosphoprotein that may upregulate the PDGFA-stimulated growth of fibroblasts and also downregulate the mitogenicity of PDGFB. The encoded protein in rodents has been shown to bind PDGFA with a low affinity.
Applications:	ELISA, IHC
Name of antibody:	PDAP1
Immunogen:	Fusion protein of human PDAP1
Full name:	PDGFA associated protein 1
Synonyms:	PAP; PAP1; HASPP28

SwissProt:	Q13442
ELISA Recommended dilution:	5000-10000
IHC positive control:	Human tonsil
IHC Recommend dilution:	50-300

