

兔抗 PDIA2 多克隆抗体

中文名称：兔抗 PDIA2 多克隆抗体

英文名称：Anti-PDIA2 rabbit polyclonal antibody

别名：PDI, PDA2, PDIP, PDIR

相关类别：一抗

储存：冷冻（-20℃）

宿主：Rabbit

抗原：PDIA2

反应种属：Human

标记物：Unconjugate

克隆类型：rabbit polyclonal

技术规格

Background:	Protein disulfide isomerases (EC 5.3.4.1), such as PDIP, are endoplasmic reticulum (ER) resident proteins that catalyze protein folding and thiol-disulfide interchange reactions
Applications:	ELISA, IHC
Name of antibody:	PDIA2
Immunogen:	Synthetic peptide of human PDIA2
Full name:	Protein disulfide isomerase family A, member 2
Synonyms :	PDI, PDA2, PDIP, PDIR
SwissProt:	Q13087
ELISA Recommended dilution:	2000-10000
IHC positive control:	Human ovarian cancer and Human lung cancer

IHC Recommend dilution:	50-200
--------------------------------	--------

